



Cubit Medisurge LLP

SURETY FOR QUALITY



PRODUCT LIST

**Manufacturer of
Disposable Medical Devices**



ABOUT COMPANY

Cubit Medisurge LLP was established as a manufacturer of disposable medical and surgical devices in 2020. We have a state-of-the-art manufacturing facility equipped with the latest technology for producing high-quality products. Our mission is to establish Cubit Medisurge as a prominent and trusted brand for medical devices and disposables. The medical equipment manufacturers in India have been critical to the country's rise as a global destination for the production, development, and export of medical equipment.

Cubit is renowned for high-quality products, and our goal is to strive towards a zero-defect approach to provide products with high precision. Microscopic care is taken at every step of the process, from sourcing raw materials to production, inspection, packaging, and sterilization.

We have a dedicated team of technical professionals who ensure that our products are designed and developed to conform to the highest international standards. We have a well-established laboratory for conducting physicochemical, microbiological, toxicity, and sterility tests. Our deep understanding of the needs of the medical fraternity comes from our other verticals of pharmaceuticals and life sciences.

Cubit Group- our journey began in 2003, and today we have more than 400 brands in India. We have expanded our Pharma product portfolio to include Cardio/ Diabetic/ CNS/ Psychiatric/ Derma/ Ophthalmology and Gynaecology divisions. We have a pan-India presence for our products with a robust marketing network of distributors, regional sales managers, marketing representatives, and C&F Agents. We export to more than 18 countries, and newer markets are being added.

CUBIT : Products known for Quality, Reliability, and Affordability

Highlights

- ISO & CE Certified Plant
- Precision Injection Molding of Components
- Modern Assembly Line and Packaging Machines
- Class 10000 and Class 100000 Clean Room
- State of the Art HVAC Systems
- Controlled environment through Hepa Filters, positive pressure, and controlled temperature
- Ethylene Oxide Sterilization [ETO]

Our Product Categories

- Blood Transfusion Therapy
- Infusion Therapy
- Gastroenterology
- Surgery
- Dialysis
- Anaesthesia
- General Purpose
- Urology

Opportunities:

- OEM /Contract Manufacturing.
- Supply of Components for Medical Devices.
- Custom-made products for Buyer/Institutions/Hospitals.

INTRAVENOUS. (I.V. CANNULA)



**I V Cannula with Wings
and with Injection port**

Size - 18G, 20G, 22G, 24G



**I V Cannula with Wings
without Injection port**

Size - 18G, 20G, 22G, 24G



**I V Cannula without Wings
and without Injection port**

Size -18G, 20G, 22G, 24G

SIZE	COLOR	CATHETER OD (mm)	CATHETER LENGTH (mm)	FLOW RATE (ml/min)
14G	Orange	2.2	45	270
16G	Grey	1.8	45	200
17G	White	1.5	45	125
18G	Green	1.3	45	90
20G	Pink	1.1	32	55
22G	Blue	0.9	25	36
24G	Yellow	0.7	19	18
26G	Violet	0.6	19	13

INTRAVENOUS INFUSION SETS - VENTED & NON VENTED



Dial Flow



Economy latex tube



Y connector



Premium



Standard



IV Set vented

Available Model :

1. Economy
2. Standard
3. Premium
4. Dial Flow

Options Available :

1. PVC / DEHP Free / PVC Free
2. Latex / Latex Free
3. Injection Site : Tube / Bulb
4. EO Sterile
5. HMHD / Ribbon / Paper Poly
6. Inner Bag / Inner Box

BLOOD TRANSFUSION THERAPY - VENTED & NON VENTED



Non Vented



Vented

Available Model :

1. Economy
2. Standard

Options Available :

1. PVC / DEHP FREE
2. 80mm/100mm Chamber
3. Tube / Bulb
4. Luer Tip / Luer Lock
5. HMHD / Ribbon
6. Inner Bag / Inner Box

INFUSION THERAPY



**Measure Volume Set
Vented or Non-Vented**

150ml, 110ml



**Pressure
Monitoring Line**

Length (cm): 10, 15, 25,
50, 100, 150, 200



Three Way Stop Cock

360° Rotation
45° click resistance
4.5 Bar Pressure (65 psi)



**Extension Line
With T Connector**

Length (cm): 10, 15,
25, 50, 100, 150, 200



**Syringe with Needle
& without Needle**

1ml, 2ml, 3ml,
5ml, 10ml, 20ml, 50ml



Extension Line

Length (cm): 10, 15,
25, 50, 100, 150, 200



**Extension Line With
Three Way Stop Cock**

Length (cm): 10, 15,
25, 50, 100, 150, 200



**Extension Line
With Dial Flow**

Length (cm): 100

GASTROENTEROLOGY



Stomach Tube

Size (FG): 8, 10, 12,
14, 16, 18, 20, 22, 24, 26



Infant Feeding Tube

Size (FG): 4, 5, 6, 7,
7.5, 8, 9, 10



Ryles Tube

Size (FG): 6, 8, 10, 12, 14,
16, 18, 20, 22, 24



Levin's Tube

Size (FG): 6, 8, 10, 12, 14,
16, 18, 20, 22, 24

UROLOGY



Foley Balloon Catheter (2way)

Size (FG) : 8, 10, 12, 14, 16, 18, 20, 22, 24



Foley Balloon Catheter (3way)

Size (FG) : 16, 18, 20, 22, 24



Nelaton Catheter / Urethral Catheter

Size (FG) : 6, 8, 10, 12, 14, 16, 18, 20, 22, 24



Urine bag with T valve

Capacity 2000 ml.
Tube length - 90 -100 cm
Bottom outlet - T - Type



Urine Collection Bag

Model : Bottom Out Let,
Top Out Let



Urine bag with Hanger and Urometer

Holds up to 2000 ml.
100-250 ml chamber
for precise urine output.

RAPID TEST SECTION



Malaria Rapid Diagnostic Test Kit

Detect malaria parasites in human body.



Ovulation (LH) Rapid Test Kit

Detects LH (Luteinizing Hormone) surge in urine.



One Step HCG Pregnancy Test Kit

To confirm pregnancy at home or clinic quickly and privately.



HIV 1/2 Rapid Test Kit

To screen for HIV infection. For early diagnosis to start timely treatment.



Typhoid Rapid Test Kit

To diagnose typhoid infection. To differentiate between current (IgM) and past(IgG) infections.



Syphilis Rapid Test Kit

To screen for syphilis infection — especially in early or latent stages.

ANESTHESIA



Oxygen Mask

Size : Adult, Child



Nasal Oxygen Cannula

Size : Adult, Child



Suction Catheter With Funnel

Size (FG): 4, 5, 6, 8, 10, 12, 14, 16, 18, 20, 22, 24



Endotracheal Tube With Cued

Size : 2, 2.5, 3, 3.5, 4, 4.5, 5, 5.5, 6, 6.5, 7, 7.5, 8, 8.5, 9



Nebulizer Mask

Size : Adult, Child



Nasal Oxygen Catheter

Size (FG): 10, 12, 14, 16, 18



Suction Catheter With Thumb Control

Size (FG): 4, 5, 6, 8, 10, 12, 14, 16, 18, 20, 22, 24



Endotracheal Tube Without Cued

Size : 2, 2.5, 3, 3.5, 4, 4.5, 5, 5.5, 6, 6.5, 7, 7.5, 8, 8.5, 9

SURGERY



Yankauer Suction Set

Length (cm): 200, 210, 250, 300



Cubit Medisurge LLP

SURETY FOR QUALITY

CERTIFICATES



ISO CERTIFICATE



CE CERTIFICATES



Contact No: - +91-9664672647, +91-9998006362

Email - info@cubitmedisurge.com

Web - www.cubitmedisurge.com

Unit 1 - Plot 37, Ozone Industrial Zone- 2,
Bhayla, Dist. Bavla, Ahmedabad- 382220.

Unit 2 - J -19, Orchid Industrial park, Bavla
Sanand Road, Bavla, Ahmedabad - 382220

