



**Cubit Medisurge LLP**

SURETY FOR QUALITY



# PRODUCT LIST

**Manufacturer of  
Disposable Medical Devices**



## ABOUT COMPANY

Cubit Medisurge LLP was established as a manufacturer of disposable medical and surgical devices in 2020. We have a state-of-the-art manufacturing facility equipped with the latest technology for producing high-quality products. Our mission is to establish Cubit Medisurge as a prominent and trusted brand for medical devices and disposables. The medical equipment manufacturers in India have been critical to the country's rise as a global destination for the production, development, and export of medical equipment.

Cubit is renowned for high-quality products, and our goal is to strive towards a zero-defect approach to provide products with high precision. Microscopic care is taken at every step of the process, from sourcing raw materials to production, inspection, packaging, and sterilization.

We have a dedicated team of technical professionals who ensure that our products are designed and developed to conform to the highest international standards. We have a well-established laboratory for conducting physicochemical, microbiological, toxicity, and sterility tests. Our deep understanding of the needs of the medical fraternity comes from our other verticals of pharmaceuticals and life sciences.

Cubit Group- our journey began in 2003, and today we have more than 400 brands in India. We have expanded our Pharma product portfolio to include Cardio/ Diabetic/ CNS/ Psychiatric/ Derma/ Ophthalmology and Gynaecology divisions. We have a pan-India presence for our products with a robust marketing network of distributors, regional sales managers, marketing representatives, and C&F Agents. We export to more than 18 countries, and newer markets are being added.

### *CUBIT : Products known for Quality, Reliability, and Affordability*

#### Highlights

- ISO 13485:2016 Certified Manufacturing Unit
- Precision Injection Molding of Components
- Modern Assembly Line and Packaging Machines
- Class 10000 and Class 100000 Clean Rooms
- State of the Art HVAC Systems
- Controlled environment through Hepa Filters, positive pressure, controlled temperature and regulated humidity
- Ethylene Oxide Sterilization [ETO]

#### Our Product Categories

- Blood Transfusion Therapy
- Infusion Therapy
- Surgery
- Anesthesia
- Urology


#### Opportunities:

- OEM /Contract Manufacturing.
- Supply of Components for Medical Devices.
- Custom-made products for Buyer/Institutions/ Hospitals.





## I.V. CANNULA WITH WINGS AND INJECTION PORT



SIZE	COLOR	CATHETER OD (mm)	CATHETER LENGTH (mm)	FLOW RATE (ml/min)
14G	Orange	2.2	45	270
16G	Grey	1.8	45	200
17G	White	1.5	45	125
18G	Green	1.3	45	90
20G	Pink	1.1	32	55
22G	Blue	0.9	25	36
24G	Yellow	0.7	19	18
26G	Violet	0.6	19	13

### Features Available:

1. Radiopaque
2. Pyrogen Free
3. Non Toxic
4. ETO Sterile
5. Latex Free
6. Blister Packing

## BLOOD TRANSFUSION THERAPY - VENTED & NON VENTED



**Non Vented**



**Vented**

### Available Model :

1. Economy
2. Premium

### Options Available :

1. PVC / DEHP FREE
2. 80mm/100mm Chamber
3. Tube / Bulb
4. Luer Tip / Luer Lock
5. HMHD / Ribbon
6. Inner Bag / Inner Box



## INTRAVENOUS INFUSION SETS - VENTED & NON VENTED



**Dial Flow**



**Economy latex tube**

Available Model :

1. Economy
2. Y Connector
3. Premium
4. Dial Flow



**Y connector**



**Premium**

Options Available :

1. PVC / DEHP Free / PVC Free
2. Latex / Latex Free
3. Injection Site : Tube / Bulb
4. EO Sterile
5. HMHD / Ribbon / Paper Poly
6. Inner Bag / Inner Box




**Standard**



**IV Set vented**

## ANESTHESIA



**Suction Catheter With Funnel**

Size (FG): 4, 5, 6, 8, 10, 12, 14, 16, 18, 20, 22, 24




**Nasal Oxygen Cannula**

Size : Adult, Child



**Endotracheal Tube With Cued**

Size : 2, 2.5, 3, 3.5, 4, 4.5, 5, 5.5, 6, 6.5, 7, 7.5, 8, 8.5, 9



**Endotracheal Tube Without Cued**

Size : 2, 2.5, 3, 3.5, 4, 4.5, 5, 5.5, 6, 6.5, 7, 7.5, 8, 8.5, 9




## INFUSION THERAPY



### Measure Volume Set Vented or Non-Vented

150ml, 110ml



### Extension Line

Length (cm): 10, 15,  
25, 50, 100, 150, 200



### Pressure Monitoring Line

Length (cm): 10, 15, 25,  
50, 100, 150, 200



### Extension Line With Three Way Stop Cock


Length (cm): 10, 15,  
25, 50, 100, 150, 200



### Three Way Stop Cock

360° Rotation  
45° click resistance  
4.5 Bar Pressure (65 psi)

## UROLOGY & SURGERY



### Foley Balloon Catheter (2way)

Size (FG) : 8, 10, 12, 14, 16,  
18, 20, 22, 24  
Latex & Silicone



### Foley Balloon Catheter (3way)

Size (FG) : 16, 18, 20,  
22, 24  
Latex & Silicone



### Urine bag with T valve

Capacity 2000 ml.  
Tube length - 90 -100 cm  
Bottom outlet - T - Type



### Yankauer Suction Set

Length (cm): 200, 210,  
250, 300

